

cDNA SYNTHESIS

Reverse Transcription (Using SuperScript™ II RT – Invitrogen, Cat. No. 18064-014)

A 20- μ L reaction volume can be used for 1 ng–5 μ g of total RNA or 1–500 ng of mRNA.

1. Add the following components to a nuclease-free microcentrifuge tube:

- 1 μ L 50–250 ng random primers
- x μ L 1 ng to 5 μ g total RNA
- 1 μ L dNTP Mix (10 mM each)

Add sterile, distilled water to bring total to 12 μ L

2. Heat mixture to 65°C for 5 min and quickly chill on ice. Collect the contents of the tube by brief centrifugation and add (note for >4 samples, prepare a master mix):

- 4 μ L 5X First-Strand Buffer
- 2 μ L 0.1 M DTT
- 1 μ L RNaseOUT™ (40 units/ μ L) (OPTIONAL)*

*RNaseOUT™ (Cat. No. 10777-019) is required if using <50 ng starting RNA.

3. Mix contents of the tube gently and incubate at 25°C for 2 min.

4. Add 1 μ L (200 units) of SuperScript™ II RT and mix by pipetting gently up and down.

5. Incubations (RT1 program on our machine):

- 25°C for 10 min.
- 42°C for 50 min.
- Inactivate the reaction by heating at 70°C for 15 min.

6. Store at -20°C.